

Antenatal Diagnostics – precise, gentle, and low-risk. Convincing details, convincing innovation.

Antenatal medicine is concerned with pregnant women and the unborn. The antenatal diagnostic test is a fundamental component of antenatal medicine. It includes all diagnostic measures by which morphological, structural, functional, chromosomal, and molecular disorders can be recognized or excluded prior to birth.

SOMATEX® offers special instruments for the various invasive procedures of antenatal diagnostics. Due to the great strains and risks inherent to invasive antenatal checkups, easy handling as well as highest material and processing quality are crucial.

Antenatal Puncture Needle



Use

The **antenatal puncture needle** is used for the various sampling techniques of

Amniocentesis

Chordocentesis

Chorionic villus sampling



Highly rigid metal alloy

Specially developed cannula tip with maximum echogenicity

Sharp tip for easy and atraumatic penetration

Silicone coated cannula shaft





Echogenic needle tip

Antenatal Puncture Needle				
Order no.	Gauge	Diameter	Length	
190 010	22	0.7 mm	90 mm	
190 012	22	0.7 mm	120 mm	
190 014	22	0.7 mm	150 mm	
190 020	20	0.9 mm	90 mm	
190 022	20	0.9 mm	120 mm	
190 024	20	0.9 mm	150 mm	
190 032	18	1.2 mm	120 mm	

Transabdominal Chorionic Sampling Set

The SOMATEX® Transabdominal
Chorionic Sampling Set is a 2-part,
coaxial puncture system for chorionic villus sampling. The cells are taken from the
abdominal wall.

Advantages at a glance

Coaxial, 2-part puncture system for optimal needle coordination

Highly rigid metal alloy

Echogenic needle tip with ultracut

Atraumatic penetration of the abdominal wall with extremely sharp cannula tip



Cannula with cm marking



Transabdominal Chorionic Sampling Set					
Order no.	Gauge	Diameter	Length		
800 155	19.5	0.95 mm	150 mm		
	22	0.7 mm	220 mm		
800 157	17.5	1.3 mm	150 mm		
	19.5	0.95 mm	220 mm		

Transcervical Chorionic Sampling Set

Advantages at a glance

Flexible catheter made of low-friction, biocompatible plastic

Good guidance through shapeable metal obturator

Sharp catheter tip

Echogenic and thereby visible in ultrasound


Sampling Set is a biopsy catheter for chorionic villus sampling. Cells are taken transcervically from the portio in the uterus.

The **SOMATEX® Transcervical Chorionic**

Catheter

Luer connector

Handle



Antenatal Diagnostics

Antenatal Diagnostics

Metal obturator

Intrauterine Drainage Set

The **SOMATEX®** Intrauterine Drainage **Set** serves for the treatment of fetuses
with obstructive urinary tract disorders.
Further indications are the treatment of
kidney congestion, bladder malfunction,
and fluid accumulation in the abdominal
cavity.

A major advantage of this product is that it can be repositioned after a wrong or not optimal placing of the shunt in the respective target area.

Advantages at a glance

Extremely sharp needle tip (trocar)

Special shunt design for secure fixation

Correctable stent

Trocar needle tip

Thread to repositioning and recovery

Puncture cannula

Lateral holes to drain fluid

Intrauterine Drainage Set

Order no. 410 055

Diameter 4 French

Obturator

ntenatal Diagnostics

Antenatal Diagnostic





Rheinstraße 7d D-14513 Teltow Germany

Tel. +49 (0) 3328 30 76-0 Fax +49 (0) 3328 30 76-99 info@somatex.com www.somatex.com