

# QUICKVUE<sup>®</sup>+

Infectious Mononucleosis TEST

## Test and Treat Today!

**Announcing the return** of QuickVue+ Infectious Mononucleosis, a rapid test for the detection of infectious mononucleosis IgM antibodies in whole blood, plasma and serum. IgM antibodies appear in early acute phase infection and disappear within 4 to 6 weeks.

QuickVue+ Infectious Mononucleosis has no age limitations and combines the simplicity of a rapid test with accurate results you can trust. The clearly distinguishable plus or minus readout and simple procedure makes it easy to run in your lab or office.



### Rapid Diagnosis and Treatment

- Results in about 5 minutes.
- Test and treat in the same office visit.
- Reduce the spread of germs.

### Easy to Use

- Self-contained Test Cassette.
- Easy-to-read plus/minus (+/-) readout.
- Use with multiple sample types.
- CLIA waived for whole blood.
- No age restrictions.

### Accurate Detection

- IgM-specific test.
- Sensitivity:  $\geq 99\%$ : Specificity:  $> 99\%$ .
- $> 99\%$  agreement with slide hemagglutination and EIA.

### Convenient Packaging

- Room temperature storage with up to 24-month shelf life from date of manufacture.
- Built-in Positive and Negative procedural controls ensure accuracy.
- Positive and Negative external controls provided.


**QUIDEL<sup>®</sup>**

For more information on QuickVue+ Infectious Mononucleosis test kits, visit our website or call your Quidel Account Manager at 1.800.874.1517 (U.S. only), or 858.552.1100 (outside U.S.).


Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • [quidel.com](http://quidel.com)




## Whole blood/hanging drop procedure:



## Whole blood/capillary tube procedure:



## Venipuncture, serum or plasma procedure:



To order: Call your Quidel Account Manager or your local distributor

QuickVue+ Infectious Mononucleosis – 20 Test Kit: Catalog #20121

**QUIDEL**<sup>®</sup>

For more information on QuickVue+ Infectious Mononucleosis test kits, visit our website or call your Quidel Account Manager at 1.800.874.1517 (U.S. only), or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • [quidel.com](http://quidel.com)

Distributor Information: