


FOR SOFT TISSUE BIOPSIES

Somacore Automatic



High quality tissue samples, intelligent handling and design are the benchmark for Somacore Automatic. Excellent usability for various clinical applications and additional smart features make Somacore Automatic the perfect solution for your soft tissue biopsy procedures.

- A large sample notch and sharp cutting cannula to achieve excellent validity of tissue samples
- Simple and intuitive one-hand operation and two easy to push trigger buttons
- Improved ultrasound visibility through echogenic markings



SOMATEX Medical Technologies GmbH
Hohenzollerndamm 150 / 151 · 14199 Berlin · Germany
Tel.: +49 (0) 30 319 82 25-00 · Fax: +49 (0) 30 319 82 25-99
info@somatex.com · www.somatex.com

SOMATEX
MEDICAL TECHNOLOGIES

	Somacore Automatic			Core Biopsy Puncture Sheath		
Color	Article No.	Diameter	Length	Article No.	Diameter	Length
▲	900340	14G	100mm	170440	13G	50mm
▲	900342	14G	160mm	170442	13G	110mm
▲	900344	14G	200mm	170444	13G	150mm
▲	900346	16G	100mm	170446	15G	50mm
▲	900348	16G	160mm	170448	15G	110mm
▲	900350	16G	200mm	170450	15G	150mm
▲	900352	18G	100mm	170452	17G	50mm
▲	900354	18G	160mm	170454	17G	110mm
▲	900356	18G	200mm	170456	17G	150mm
▲	900358	18G	250mm	170458	17G	200mm
▲	900360	20G	100mm	170460	19G	50mm
▲	900362	20G	160mm	170462	19G	110mm
▲	900364	20G	200mm	170464	19G	150mm

Somacore Automatic comes in four colors, which correspond to different needle diameters.

Somacore Automatic is used for taking samples of soft tissues such as breast, liver, kidney, prostate, spleen, lymph nodes and various soft tissue tumors. It is optional available with core biopsy puncture sheaths.

