Tumark®Professional

The new marker system



SOMATEX MEDICAL TECHNOLOGIES GMBH

www.somatev.com

Tumark® Professional -

Precision visibility

Safe and easy thanks to single-handed operation

Tumark® Professional enables precise tissue marking. With Tumark® Professional, tumors can be marked in a gentle way and with minimum invasion before cytostatic or radiation treatment begins. Following successful chemo or radiation therapy, the tissue area concerned can be reliably and quickly found again. Tumark® Professional is ideally suited for marking a biopsy site.

The tumor can be reliably localised and precisely radiated - thereby protecting healthy areas. Thanks to its 3D design, the marker offers firm anchorage in the tissue. Moreover the marker's surface characterisitcs ensure optimal visibility for ultrasound, MRI and mammography.

The Tumark® Professional is also available as fully MRI-compatible version: Tumark® MRI.

Advantages at a glance

- Clever 3D marker design enables firm anchorage in the tissue and increases visibility in all positions
- Optimal visibility for ultrasound, mammography and MRI
- Marker is MRI compatible up to 3 Tesla
- The marker is approved as implantable material (Nitinol)
- Ergonomic handle for single-handed operation

Use **Advantages Product options**

Use

- Marking after removal of tumor
- Marking of suspicious lesions
- Tumor marking during neoadjuvant chemotherapy
- For marking the biopsy site
- Orientation for planning radiation therapy

The marker is supplied preloaded



Ultrasound marking

Ultrasound cut



Mammography





MRI

Ultrasound





Product options				
Product name	Order-No.	Cannula diameter	Packaging unit	
Tumark® Professional	271 560	18G 1,2 x 120 mm	10	
NEW: Tumark® MRI*	601 570	18G 1,2 x 120 mm	10	

^{*} Shape and structure of the MRI version differ from figures

Images from Dr Sigamany	Padiology Diganoctics	University Hospital Tühingen
TITICALES ITOTILI DE SIEGITIANT	KAMOOMA I JAMOOMA	University mospital Hiblingen

SOMATEX® MEDICAL TECHNOLOGIES GMBH

Rheinstr. 7d D-14513 Teltow Germany

Tel. +49 (0) 3328 30 76-0 Fax +49 (0) 3328 30 76-99 info@somatex.com www.somatex.com